


PATIENT GROUPS AND APPROPRIATE FLEXINEB E2 MASK SIZE


Based on the muzzle circumference dimension the correct size Mask is chosen as follows:



Small/Foal/Pony Mask
Muzzle circumference measurement 12" to 17" can be used with foals and small ponies

Standard/Medium Mask
Muzzle circumference measurement 17" to 20" can be used with Thoroughbreds, Standardbreds, Quarter horses and small Warmbloods

Large Mask
Muzzle circumference measurement 20" to 26" can be used with large Warmbloods.



CONVERTING INCHES TO CENTIMETRES (cm)

Foal:	12 inches to 17 inches or 30.5cm to 43cm
Standard/Medium Mask (Arabian horse):	17 inches to 20 inches or 43cm to 51cm
Large Mask (Warmblood):	20 inches to 26 inches or 51cm to 66cm

Please do not hesitate to contact Catherine Trindall (0417415355), at Trindall Equestrian Services Pty Ltd if you have questions regarding the appropriate size mask to purchase.